

LiNbO₃ Crystal

WISOPTIC Capabilities of LiNbO₃

- Various sizes of finished components for different applications.
- Strict quality control
- Very competitive price and reliable delivery
- Technical support at any time


WISOPTIC Standard Specifications of LiNbO₃

| | |
|----------------------------------|--|
| Dimensions | (4×4~9×9)×25 mm |
| Dimension Tolerance | ± 0.1 mm |
| Angle Tolerance | ± 0.5° |
| Flatness | < λ/8 @ 632.8 nm |
| Surface Quality | < 20/10 [S/D] |
| Parallelism | < 20" |
| Perpendicularity | ≤ 5' |
| Chamfer | ≤ 0.2mm @ 45° |
| Transmitted Wavefront Distortion | < λ/4 @ 632.8 nm |
| Capacitance | 20 pF |
| Electrodes | Gold |
| λ/4 Voltage | ~1.5 kV (@ 1064 nm) |
| Dynamic Extinction Ratio | 20 dB (@ 1064 nm) |
| Clear Aperture | 8.5 mm (> 90% central area) |
| Coating | AR coating: R < 0.2% @ 1064 nm, R < 0.5% @ 532 nm |
| Laser Damage Threshold | LN 100 MW/cm ² @ 10ns, 10Hz, 1064nm MgO : LN 150 MW/cm ² @ 10ns, 10Hz, 1064nm |

Please note: PPLN/PPMgLN and SLN/MgSLN (Stoichiometric LN) is available.

PP-LiNbO₃ & PP-MgO:LiNbO₃


Chip Lay-out



Main Specifications

| Items | Properties | Inspection |
|--------------------------------------|--|-------------------|
| Material | 5 mol.% MgO:LN | NA |
| Period (Λ , μm) | 18.20, 18.50, 18.80, 19.10, 19.40, 19.70, 20.00, 20.30, 20.60, 20.90 | Microscope |
| Main Function | Second Harmonic Generation | NA |
| Parallelism/Perpendicularity | $\leq 5' / 21'$ | Autocollimator |
| Flatness | $\leq \lambda / 6$ ($\lambda = 633 \text{ nm}$) | Interferometer |
| Scratch/Dig | $\leq 20 / 10$ | Microscope |
| Optical coating (S1/S2 facets) | S1/S2 @ 750-810 (R < 0.5%) / 1500-1620 (R < 0.5%) nm | Spectral Analyzer |
| Aperture Size | 12.3 x 1.0 mm ² (W x T) | Cutting Machine |
| Available Length | 10/25/50 \pm 0.2 mm | |
| Channel Clear Aperture | $\geq 80\%$ (T), $\geq 90\%$ (W) | NA |

Phase Matching Tuning Curve




Oven & Temperature Controller



Oven and temperature controller


The Temperature Controller is designed to offer optimum temperature control of the PPLN crystal within the oven.

| Specifications | |
|--|---|
| Temperature Controller | Oven |
| <ul style="list-style-type: none"> • Purpose ASIC processor • PID control • RS-232 interface (optional) • Programmable (optional) • Stability: 0.1 C/±0.01 C • Size: 70×145×150mm³ • Normal Package Standard • Voltage: 220V AC | <ul style="list-style-type: none"> • Design: Radioactive heat compensation • Plane temperature distribution • Size: 38×42×90mm³ 38×54×112mm³ 33×70×110mm³ • Sensor: Pt100 thermocouple • Working temperature: ≤100 C/≤200 C |



Applications:

Heating nonlinear crystals is usually employed in NCPM, OPO, OPA, etc.



Note:

- 1、 Please confirm it when order and check it before plug in the power. Burned and other damages, which caused by improper power selection, are not guaranteed to repair.
- 2、 Special oven size, holder of crystals and right-angle support setting applications are available upon request.

Operation description:

1. Connect "AC220V" power line and "6 PIN with D89" splice to the temperature control furnace.
2. The electric power will come when you push "switch" to "—", after that the screen will display "N: 25.0°C" and "S: 25.0°C ", the "N" and "S" represent the present temperature and the setting temperature respectively.
3. Press the "set" button once, the cursor will be on hundreds of temperature numbers, you can set this number, with "Add" button you can increase the setting temperature, with "Sub" button you can decrease it.
4. Press the "set" button twice, you can set the tens of temperature numbers.
5. Press the "set" three times, you can set the single digit of temperature numbers.
6. Press the "set" four times, the temperature control furnace start to work.