



Dye Laser Handpiece 585nm

Working Mode: input 532nm,
output 585nm
Conversion Efficiency: 65% ~ 75%
Beam Diameter: input 8mm,
output 2~3mm
Operating Life: > 40k shots
Outside Dimensions: $\phi 31.5 \times 178$ mm
Size of Connector: M22 \times 1.5mm
Main Application: Treatment of facial
telangiectasia

Dye Laser Handpiece 650nm

Working Mode: input 532nm,
output 650nm
Conversion Efficiency: 45% ~ 55%
Beam Diameter: input 8mm,
output 2~3mm
Operating Life: > 40k shots
Outside Dimensions: $\phi 31.5 \times 178$ mm
Size of Connector: M22 \times 1.5mm
Main Application: Removal of green
tattoo